

Technical Data Sheet

Product	Anti-Hu CD300e Purified Low Endotoxin	
Cat. Number/Size	12-728-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD300e
Clone	UP-H2
Format	Purified Low Endotoxin
Reactivity	Human
Application	FC (QC tested), IP, WB, FUNC
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody UP-H2 recognizes an extracellular epitope on CD300e / IREM-2, a 32 kDa glycoprotein expressed by mature monocytes and peripheral blood myeloid dendritic cells.
Other names	CLM2, CMRF35-A5, LMIR6, IREM-2, PlgR2
Immunogen	CD300e-HA-transfected cells
Entrez Gene ID	342510
Gene name	CD300E
NCBI Full Gene Name	CD300e molecule
UniProt ID	Q496F6

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Low endotoxin azide free phosphate buffered saline (PBS) solution, 0.2 µm filter sterilized. Endotoxin level is less than 0.01 EU/µg of the protein, as determined by the LAL test.
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References	www.exbio.cz
-----------------------	--

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.