

Technical Data Sheet

| | | |
|------------------|---|-----------------|
| Product | Anti-Hu CD178 Purified | |
| Cat. Number/Size | 11-754-C025 | 0.025 mg |
| | 11-754-C100 | 0.1 mg |
| | For Research Use Only. | |
| | Not for use in diagnostic or therapeutic procedures. | |

| | |
|---------------------|---|
| Antigen | CD178 |
| Clone | NOK-1 |
| Format | Purified |
| Reactivity | Human |
| Application | FC (QC tested), IP, WB, ICC |
| Application details | Flow cytometry: Recommended dilution: 1-4 µg/ml |
| Isotype | Mouse IgG1 |
| Specificity | The mouse monoclonal antibody NOK-1 recognizes an extracellular epitope of CD178 / Fas-L, an approximately 40 kDa transmembrane glycoprotein expressed on neutrophils, monocytes, and activated T and NK cells. |
| Other names | APTL, FASL, FAS-ligand, CD95L, TNFSF6, TNLG1A |
| Workshop | HLDA VII: 70322 |
| Immunogen | L5178Y mouse T lymphoma cells expressing recombinant human CD178 |
| Entrez Gene ID | 356 |
| Gene name | FASLG |
| NCBI Full Gene Name | Fas ligand |
| UniProt ID | P48023 |

| | |
|----------------------|---|
| Concentration | 1 mg/ml |
| Preparation | Purified by protein-A affinity chromatography. |
| Formulation | Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide |
| Storage and handling | Store at 2-8°C. Do not freeze. |

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.