



	Technical Data Sheet	
Product	Anti-TIAR Purified	
Cat. Number/Size	11-920-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	TIAR	
Clone	6E3	
Format	Purified	
Reactivity	Human, Mouse, Non-human primates	
Application	WB, IHC-P, ICC, FC (QC tested), ELISA	
Application details	Flow cytometry: Intracellular staining; recommended dilution: 5 µg/ml	
Isotype	Mouse IgG2a kappa	
Specificity	The mouse monoclonal antibody 6E3 recognizes TIAR, an ubiquitously expressed approximately 42 kDa intracellular RNA-binding protein.	
Other names	TCBP, TIAL1	
Immunogen	GST-tagged recombinant human TIAR	
Entrez Gene ID	7073	
Gene name	TIAL1	
NCBI Full Gene Name	TIA1 cytotoxic granule associated RNA binding prot	
UniProt ID	Q01085	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.