

Technical Data Sheet

Product	Anti-Hu CD61 Purified	
Cat. Number/Size	11-647-C025	0.025 mg
	11-647-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD61
Clone	VIPL2
Format	Purified
Reactivity	Human, Non-human primates
Application	FC (QC tested), WB, IHC(F)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody VIPL2 recognizes an extracellular epitope of CD61, a 90-110 kDa transmembrane glycoprotein of integrin family, expressed on platelets, megacaryocytes, osteoclasts, endothelial cells and other cell types, including leucocytes and smooth muscle cells.
Other names	integrin beta 3, ITGB3, GP3A
Workshop	HLDA V: WS Code 5T-124
Entrez Gene ID	3690
Gene name	ITGB3
NCBI Full Gene Name	integrin subunit beta 3
UniProt ID	P05106

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.