



	Technical Data Sheet	
Product	Anti-Hu CD318 Purified	
Cat. Number/Size	11-833-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD318	
Clone	CUB1	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 2-6 µg/ml	
Isotype	Mouse IgG2b	
Specificity	The mouse monoclonal antibody CUB1 recognizes an extracellular epitope of CD318, a type I transmembrane protein involved in early hematopoiesis.	
Other names	CDCP1, SIMA135, TRASK	
Workshop	HCDM listed	
Immunogen	NIH-3T3/CD318 cells	
Entrez Gene ID	64866	
Gene name	CDCP1	
NCBI Full Gene Name	CUB domain containing protein 1	
UniProt ID	Q9H5V8	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.