

Technical Data Sheet

Product	<b>Anti-AGR3 Purified</b>	
Cat. Number/Size	<b>11-560-C025</b>	<b>0.025 mg</b>
	<b>11-560-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

Antigen	AGR3
Clone	AGR3.1
Format	Purified
Reactivity	Human
Application	WB, IHC(P), IHC(F), ICC
Isotype	Mouse IgG1
Specificity	The antibody AGR3.1 recognizes the epitope HETTDKNLS within the AGR3 (AG3) protein (19-20 kDa); a secreted cytoplasmic protein which can serve as a marker of carcinogenesis.
Other names	AG3, HAG3
Immunogen	Purified human AGR3 protein
Entrez Gene ID	155465
Gene name	AGR3
NCBI Full Gene Name	anterior gradient 3, protein disulphide isomerase
UniProt ID	Q8TD06

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References [www.exbio.cz](http://www.exbio.cz)

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.