



	Technical Data Sheet	
Product	Anti-ESAT-6 (M. tuberculosis) Purified	
Cat. Number/Size	11-491-C025	0.025 mg
	11-491-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	ESAT-6 (M. tuberculosis)	
Clone	PAb (491)	
Format	Purified	
Reactivity	Mycobacterium tuberculosis	
Application	WB	
Application details	Western blotting: Recommended dilution 0.5-1 µg/ml, reducing conditions.	
Isotype	Rabbit polyclonal	
Specificity	The rabbit polyclonal antibody to EsaT-6 recognizes the EsaT-6 (Rv3875) protein of Mycobacterium tuberculosis and Mycobacterium bovis.	
Other names	ESXA, ESAT-6, Rv3875	
Immunogen	EsaT-6 (Rv3875) of <i>Mycobacterium tuberculosis</i>	
Entrez Gene ID	886209	
Gene name	esxA	
NCBI Full Gene Name	ESAT-6 protein EsxA	
UniProt ID	P9WVK7	
Concentration	1 mg/ml	
Preparation	Purified by ligand affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.