

Technical Data Sheet

| | | |
|------------------|-------------------------------------------------------------|---------------|
| Product | Anti-Blood Group B Culture supernatant | |
| Cat. Number/Size | 99-806-L001 | 1.0 ml |
| | For Research Use Only. | |
| | Not for use in diagnostic or therapeutic procedures. | |

| | |
|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Antigen | Blood Group B |
| Clone | HEB-29 |
| Format | Culture supernatant |
| Reactivity | Human |
| Application | AGG, IHC(P) |
| Isotype | Mouse IgM |
| Specificity | The antibody HEB-29 reacts with human blood group B. The specificity of the antibody HEB-29 was confirmed by comparison of specificity and reactivity to standard reagent using >5000 samples of blood. |
| Immunogen | Mixture of erythrocytes of group B and glycoprotein fraction isolated from saliva of secretors with blood group B. |

| | |
|----------------------|-------------------------------------------------------------------------|
| Preparation | Hybridoma culture supernatant concentrated by membrane ultrafiltration. |
| Storage and handling | Store at 2-8°C. Do not freeze. |

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.