



	Technical Data Sheet	
Product	Anti-Hu CD184 APC	
Cat. Number/Size	1A-146-T100	100 tests
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD184	
Clone	12G5	
Format	APC	
Reactivity	Human, Non-human primates	
Application	FC (QC tested)	
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.	
Excitation laser	red (633 nm)	
Isotype	Mouse IgG2a kappa	
Specificity	The mouse monoclonal antibody 12G5 recognizes an extracellular epitope of CD184, a 45 kDa G-protein-linked CXC chemokine receptor widely expressed on blood and tissue cells.	
Other names	CXCR4, LAP3, LCR1, HM89, fusin, SDF1R	
Workshop	HLDA VII: WS Code 70204	
Immunogen	CP-MAC-infected Sup-T1 cells	
Entrez Gene ID	7852	
Gene name	CXCR4	
NCBI Full Gene Name	C-X-C motif chemokine receptor 4	
UniProt ID	P61073	
Preparation	Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.