

Technical Data Sheet

Product

Anti-Hu CD270 APC

Cat. Number/Size

1A-911-T100

100 tests

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CD270
Clone	CW10
Format	APC
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Excitation laser	red (633 nm)
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody CW10 recognizes an extracellular epitope on CD270, a type I transmembrane protein expressed on resting T cells, monocytes, and immature dendritic cells.
Other names	TNFRSF14, TR2, ATAR, HVEA, HVEM, LIGHTR
Immunogen	recombinant human CD270
Entrez Gene ID	8764
Gene name	TNFRSF14
NCBI Full Gene Name	TNF receptor superfamily member 14
UniProt ID	Q92956

Preparation Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Formulation Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage and handling Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.