

Technical Data Sheet

Product	Anti-Hu CD141 Purified	
Cat. Number/Size	11-795-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD141
Clone	M80
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody M80 recognizes an extracellular epitope of CD141, a 75 kDa transmembrane glycoprotein expressed mainly on macrophages, monocytes, platelets and endothelial cells.
Other names	Thrombomodulin, fetomodulin, THBD, THRM, THPH12, TM
Immunogen	MV4-11 cell line
Entrez Gene ID	7056
Gene name	THBD
NCBI Full Gene Name	thrombomodulin
UniProt ID	P07204

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.