

Technical Data Sheet

Product	<b>Anti-Hu CD5 Purified</b>	
Cat. Number/Size	<b>11-691-C025</b>	<b>0.025 mg</b>
	<b>11-691-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	
Antigen	CD5	
Clone	L17F12	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP, WB, ICC, IHC	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml. Western blotting: Laurylmaltoside lysing buffer; non-reducing conditions; recommended dilution: 1-2 µg/ml.	
Isotype	Mouse IgG2a kappa	
Specificity	The mouse monoclonal antibody L17F12 reacts with an extracellular epitope of CD5, a 67kDa single-chain transmembrane glycoprotein expressed on mature T lymphocytes, most of thymocytes and B lymphocytes subset (B-1a lymphocytes).	
Other names	T1, LEU1	
Immunogen	Human acute lymphoblastic leukemia (ALL) T cells	
Entrez Gene ID	921	
Gene name	CD5	
NCBI Full Gene Name	CD5 molecule	
UniProt ID	P06127	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	<a href="http://www.exbio.cz">www.exbio.cz</a>	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.