



Technical Data Sheet

Product	Anti-Hu p53 (pS392) Purified	
Cat. Number/Size	11-467-C025	0.025 mg
	11-467-C100	0.1 mg
For Research Use Only.		
Not for use in diagnostic or therapeutic procedures.		
Antigen	p53 (pS392)	
Clone	FP3.2 [FPS392]	
Format	Purified	
Reactivity	Human	
Application	WB, IHC-P	
Application details	Immunohistochemistry (paraffin sections): Standard ABC technique (DAB+), pretreatment: high temperature antigen retrieval (microwave, pressure cooker) in 10 mM citrate buffer pH 6.0 or 1 mM EDTA-NaOH buffer pH 8.0, recommended dilution: 10 µg/ml, incubation: 1 hour at RT; or overnight at 4°C, positive tissue: breast carcinoma with high level of wild-type p53. Western blotting: recommended dilution: 1 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody FP3.2 [FPS392] reacts with human p53 tumour suppressor intracellular protein phosphorylated at CKII site (Ser 392).	
Other names	BCC7, TRP53, TP53, LFS1	
Immunogen	KLH-conjugated phosphopeptide RHKKLMFKTEGPDS[P]D, corresponding to amino acids 378-393 of human p53.	
Entrez Gene ID	7157	
Gene name	TP53	
NCBI Full Gene Name	tumor protein p53	
UniProt ID	P04637	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.