

Technical Data Sheet

Product	Anti-Biotin Purified
Cat. Number/Size	11-944-C100 0.1 mg
	For Research Use Only.
	Not for use in diagnostic or therapeutic procedures.

Antigen	Biotin
Clone	1D4-C5
Format	Purified
Application	FC (QC tested), ELISA, IP, MC (CyTOF)
Isotype	Mouse IgG2a kappa
Specificity	Mouse monoclonal antibody 1D4-C5 recognizes biotin and can be used as a second-step reagent for detecting biotinylated antigens.
Other names	vitamin H, vitamin B7
Immunogen	biotinylated limulus hemocyanin

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.