



	Technical Data Sheet	
Product	Anti-CRP Purified	
Cat. Number/Size	11-537-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CRP	
Clone	C2	
Format	Purified	
Reactivity	Human	
Application	ELISA, IHC-P	
Application details	ELISA: The antibody C2 has been tested as the capture antibody in a sandwich ELISA for analysis of human CRP.	
Isotype	Mouse IgG1	
Specificity	The mouse monoclonal antibody C2 recognizes human CRP, an 117 kDa plasma acute phase protein; it recognizes the antigen both in the presence and in the absence (samples containing EDTA) of Ca ²⁺ .	
Immunogen	Purified human CRP	
Entrez Gene ID	1401	
Gene name	CRP	
NCBI Full Gene Name	C-reactive protein	
UniProt ID	P02741	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.