

Technical Data Sheet

Product	Anti-Hu CD64 Biotin	
Cat. Number/Size	1B-644-C025	0.025 mg
	1B-644-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD64
Clone	10.1
Format	Biotin
Reactivity	Human, Non-human primates
Application	FC (QC tested), IP, WB, IHC(P), IHC(F), ICC
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody 10.1 recognizes an extracellular epitope on CD64/FcgammaRI, a 72 kDa single chain type I glycoprotein, that is expressed on monocytes/macrophages, dendritic cells, and activated granulocytes.
Other names	FcRI, IGFR1, FcγR1A
Workshop	HLDA III: WS Code M-250
Immunogen	Rheumatoid synovial fluid cells and fibronectin purified human monocytes
Entrez Gene ID	2209
Gene name	FCGR1A
NCBI Full Gene Name	Fc fragment of IgG receptor Ia
UniProt ID	P12314

Concentration	1 mg/ml
Preparation	The purified antibody is conjugated with biotin-LC-NHS under optimum conditions. The reagent is free of unconjugated biotin.
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.