

Technical Data Sheet

Product	Anti-PSMA Purified	
Cat. Number/Size	11-539-C025	0.025 mg
	11-539-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	PSMA
Clone	GCP-05
Format	Purified
Reactivity	Human
Application	FC (QC tested), IP, ICC, MC (CyTOF)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody GCP-05 recognizes extracellular domain (preferentially in native form) of glutamate carboxypeptidase II (NAALADase, FOLH1, PSMA), an approximately 95-110 kDa transmembrane glycoprotein expressed mainly in tumour neovasculatures, nervous system and jejunum, which is an important prostate tumour marker.
Other names	GCP2, FOLH1, NAALADase I, PGGCP, FGGCP, FGCP,
Immunogen	amino acids 44-750 of human GCPII
Entrez Gene ID	2346
Gene name	LOLH1
NCBI Full Gene Name	folate hydrolase 1
UniProt ID	Q04609

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.