

Technical Data Sheet

Product **Anti-Hu Tomm34 Purified**
Cat. Number/Size **11-891-C100** **0.1 mg**
For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

| | |
|---------------------|--|
| Antigen | Tomm34 |
| Clone | Tomm34-4.1 |
| Format | Purified |
| Reactivity | Human |
| Application | WB, ICC, IP |
| Isotype | Mouse IgG1 |
| Specificity | The mouse monoclonal antibody Tomm34-4.1 recognizes an epitope within TPR2 domain of human Tomm34 protein. |
| Other names | URCC3, HTOM34P |
| Immunogen | Recombinant human Tomm34 |
| Entrez Gene ID | 10953 |
| Gene name | TOMM34 |
| NCBI Full Gene Name | Translocase of outer mitochondrial membrane 34 |
| UniProt ID | Q15785 |

| | |
|----------------------|---|
| Concentration | 1 mg/ml |
| Preparation | Purified by protein-A affinity chromatography. |
| Formulation | Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide |
| Storage and handling | Store at 2-8°C. Do not freeze. |

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.