



	Technical Data Sheet	
Product	Anti-Ms CD16/CD32 PE	
Cat. Number/Size	1P-141-C025	0.025 mg
	1P-141-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD16/CD32	
Clone	93	
Format	PE	
Reactivity	Mouse	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 2 µg/ml.	
Excitation laser	blue (488 nm)	
Isotype	Rat IgG2a lambda	
Specificity	The rat monoclonal antibody 93 recognizes a common extracellular epitope of murine CD16 (FcgammaRIII) and CD32 (FcgammaRII), the low affinity receptors for IgG.	
Other names	CD16, CD32, FcgammaRIII, FcgammaRII, FCGR3, FCGR2	
Immunogen	Murine pre-B cells	
Concentration	0.5 mg/ml	
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.