

Technical Data Sheet

Product **Anti-Hu CD355 Purified**
Cat. Number/Size **11-910-C100** **0.1 mg**
For Research Use Only.
Not for use in diagnostic or therapeutic procedures.

Antigen	CD355
Clone	Cr24.1
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 3-6 µg/ml.
Isotype	Mouse IgG2b
Specificity	The mouse monoclonal antibody Cr24.1 recognizes an extracellular epitope of human CD355, a type I transmembrane glycoprotein, expressed on subsets of CD8+ and CD4+ T cells, on activated NK cells and NKT cells.
Other names	CRTAM
Immunogen	CD355-P815 cells
Entrez Gene ID	56253
Gene name	CRTAM
NCBI Full Gene Name	cytotoxic and regulatory T cell molecule
UniProt ID	O95727

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.