



	Technical Data Sheet	
Product	Anti-Hu CD274 Purified	
Cat. Number/Size	11-177-C025	0.025 mg
	11-177-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD274	
Clone	29E.2A3	
Format	Purified	
Reactivity	Non-human primates, Human	
Application	FC (QC tested), IHC-F	
Application details	Immunohistochemistry (frozen sections): Acetone fixation. Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Mouse IgG2b kappa	
Specificity	The mouse monoclonal antibody 29E.2A3 recognizes an extracellular epitope of CD274 / PD-L1 (also known as B7-H1), a 40 kDa type I transmembrane protein expressed by dendritic cells, activated T cells, activated monocytes, and in various tissues, above all in heart and skeletal muscle, placenta and lung, and in many cancer cells, including T cell lymphomas, melanomas, and glioblastomas.	
Other names	B7H1, PDL1, PDCD1L1, PDCD1LG1, PDCD1 ligand 1, B7-H1	
Immunogen	Full length human CD274	
Entrez Gene ID	29126	
Gene name	CD274	
NCBI Full Gene Name	CD274 molecule	
UniProt ID	Q9NZQ7	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.