

Technical Data Sheet

Product	Anti-TFG Purified	
Cat. Number/Size	11-709-C025	0.025 mg
	11-709-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	TFG
Clone	TFG-03
Format	Purified
Reactivity	Human
Application	IP, WB, ICC
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG2b
Specificity	The mouse monoclonal antibody TFG-03 recognizes TFG, an approximately 50 kDa intracellular protein with regulatory functions.
Other names	TF6, TRKT3, SPG57
Immunogen	Human recombinant protein Trx-his-NTFG
Entrez Gene ID	10342
Gene name	HMSNP
NCBI Full Gene Name	trafficking from ER to golgi regulator
UniProt ID	Q92734

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.