

Technical Data Sheet

Product

Anti-Hu IgG (Fab) Purified

Cat. Number/Size

11-319-C025 **0.025 mg**
11-319-C100 **0.1 mg**

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	IgG (Fab)
Clone	4A11
Format	Purified
Reactivity	Human
Application	FC (QC tested), WB, ELISA
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG2b
Specificity	The mouse monoclonal antibody 4A11 reacts with Fab fragment of human IgG. This antibody detects Fab fragments regardless light chain type.
Other names	Immunoglobulin G Fab fragment
Immunogen	Purified human IgG.

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.