

Technical Data Sheet

Product	Anti-HSP90 beta Purified	
Cat. Number/Size	11-534-C025	0.025 mg
	11-534-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	HSP90 beta
Clone	MBH90B
Format	Purified
Reactivity	Cow, Human, Mouse
Application	IP, WB, ICC
Application details	Western blotting: Recommended dilution: 1 µg/ml. Immunocytochemistry: Recommended dilution: 1 µg/ml. Immunoprecipitation: Recommended dilution: 1 µg/ml.
Isotype	Mouse IgG1
Specificity	The antibody MBH90B recognizes the EEVHHG epitope within the N-terminal part of Hsp90 beta an ubiquitously expressed intracellular protein with calculated Mw of 83.3 kDa, however, migrating as a 90 kDa band under reducing SDS-PAGE conditions.
Other names	HSP90B, HSPCB
Immunogen	Peptide corresponding to the EEVHHGEEEEVEC sequence within N-terminal part of human Hsp90.
Gene name	HSP90 beta
NCBI Full Gene Name	heat shock protein 90 beta
UniProt ID	P08238

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.