

Technical Data Sheet

Product	<b>Anti-Hu CD264 FITC</b>	
Cat. Number/Size	<b>1F-519-C025</b>	<b>0.025 mg</b>
	<b>1F-519-C100</b>	<b>0.1 mg</b>
	<b>For Research Use Only.</b>	
	<b>Not for use in diagnostic or therapeutic procedures.</b>	

---

Antigen	CD264
Clone	TRAIL-R4-01
Format	FITC
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1
Specificity	The antibody TRAIL-R4-01 reacts with an extracellular epitope of TRAIL-R4, a 42 kDa transmembrane protein expressed on various blood cells.
Other names	TNFRSF10D, DcR2, TRUNDD, TRAILR4
Workshop	HLDA IX
Immunogen	TRAIL-R4 (aa 1-210) - hlgGhc fusion protein
Entrez Gene ID	8793
Gene name	TNFRSF10D
NCBI Full Gene Name	TNF receptor superfamily member 10d
UniProt ID	Q9UBN6

---

Concentration	0.1 mg/ml
Preparation	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

---

Images and References	<a href="http://www.exbio.cz">www.exbio.cz</a>
-----------------------	--

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.