

Technical Data Sheet

Product **Anti-Hu CD150 Purified Azide Free**

Cat. Number/Size **10-660-C025** **0.025 mg**
10-660-C100 **0.1 mg**

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CD150
Clone	SLAM.4
Format	Purified Azide Free
Reactivity	Human
Application	FC (QC tested), IP, ICC
Application details	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody SLAM.4 recognizes an extracellular epitope of CD150, a cell surface molecule expressed on lymphocytes and involved in their activation.
Other names	SLAM, SLAMF1
Immunogen	Human CD150-transfected 300.19 cells
Entrez Gene ID	6504
Gene name	SLAMF1
NCBI Full Gene Name	signaling lymphocytic activation molecule family m
UniProt ID	Q13291

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.