



	Technical Data Sheet	
Product	Anti-Hu CD334 Purified	
Cat. Number/Size	11-591-C025	0.025 mg
	11-591-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD334	
Clone	4FR6D3	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), ICC	
Application details	Flow cytometry: Recommended dilution: 2-6 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody 4FR6D3 reacts with an extracellular epitope of CD334, the fibroblast growth factor receptor 4, which is an approximately 88 kDa receptor tyrosine kinase expressed in variety of tissues.	
Other names	FGFR4, TKF, JTK2, FGFR-4	
Workshop	HCDM listed	
Immunogen	NIH 3T3 cells transfected with full length human CD334	
Entrez Gene ID	2264	
Gene name	FGFR4	
NCBI Full Gene Name	fibroblast growth factor receptor 4	
UniProt ID	P22455	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.