



	Technical Data Sheet	
Product	Anti-GlcB (M. tuberculosis) Purified	
Cat. Number/Size	11-192-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	GlcB (M. tuberculosis)	
Clone	PAb (192)	
Format	Purified	
Reactivity	Mycobacterium tuberculosis	
Application	WB, ELISA	
Isotype	Rabbit polyclonal	
Specificity	The rabbit polyclonal antibody to GlcB reacts with Mycobacterium tuberculosis protein GlcB (malate synthase, Rv1837c).	
Other names	Rv1837c	
Immunogen	Recombinant GlcB protein (produced in E. coli)	
Entrez Gene ID	885713	
Gene name	glcB	
NCBI Full Gene Name	malate synthase	
UniProt ID	P9WK17	
Concentration	1 mg/ml	
Preparation	Purified by ligand affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.