

Technical Data Sheet

Product

Anti-HLA-B7 Purified

Cat. Number/Size

11-547-C025

0.025 mg

11-547-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen

HLA-B7

Clone

BB7.1

Format

Purified

Reactivity

Non-human primates, Human

Application

FC (QC tested)

Application details

Flow cytometry: Recommended dilution: 1-4 µg/ml.

Isotype

Mouse IgG1

Specificity

The mouse monoclonal antibody BB7.1 recognizes an extracellular antigen of HLA-B7 antigen. Although highly specific, it can cross-react with HLA-B42 antigen.

Immunogen

Papain solubilised HLA-A2, B7

Concentration

1 mg/ml

Preparation

Purified by protein-A affinity chromatography.

Formulation

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage and handling

Store at 2-8°C. Do not freeze.

Images and References

www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.