



	Technical Data Sheet	
Product	Anti-Hu CD209 Purified	
Cat. Number/Size	11-640-C025	0.025 mg
	11-640-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD209	
Clone	UW60.1	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IHC, ICC	
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.	
Isotype	Mouse IgG2a	
Specificity	The mouse monoclonal antibody UW60.1 recognizes an extracellular epitope of human CD209 (DC-SIGN), a 44 kDa transmembrane receptor, expressed on the surface of dendritic cells and macrophages.	
Other names	CLEC4L, DC-SIGN1, CDSIGN	
Immunogen	CD209-His-tagged fusion protein	
Entrez Gene ID	30835	
Gene name	CD209	
NCBI Full Gene Name	CD209 molecule	
UniProt ID	Q9NNX6	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.