

Technical Data Sheet

Product	Anti-Hu CD307a Purified
Cat. Number/Size	11-860-C100 0.1 mg
	For Research Use Only.
	Not for use in diagnostic or therapeutic procedures.

Antigen	CD307a
Clone	E3
Format	Purified
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody E3 recognizes an extracellular epitope of CD307a, a transmembrane glycoprotein expressed predominantly on the surface of mature B cells.
Other names	FcRL1, FcRH1, IRTA5
Workshop	HLDA IX
Immunogen	DNA-immunization followed by a boost with the CD307a transfected cells
Entrez Gene ID	115350
Gene name	FCRL1
NCBI Full Gene Name	Fc receptor like 1
UniProt ID	Q96LA6

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.