



	Technical Data Sheet	
Product	Anti-Ms CD41 Purified Low Endotoxin	
Cat. Number/Size	12-763-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD41	
Clone	MWReg30	
Format	Purified Low Endotoxin	
Reactivity	Mouse	
Application	FC (QC tested), IP, IHC-P, IHC-F, ELISA, FUNC	
Application details	Functional application: Platelet depletion. Flow cytometry: Recommended dilution: 1-4 µg/ml	
Isotype	Rat IgG1 kappa	
Specificity	The rat monoclonal antibody MWReg30 recognizes an extracellular epitope of CD41 (GPIIb), a transmembrane glycoprotein (integrin family) composed of two chains GPIIb alpha (heavy chain; 120 kDa) and GPIIb beta (light chain; 23 kDa). CD41 is mainly expressed on platelets and megakaryocytes.	
Other names	Platelet GPIIb, Integrin alpha-IIb, GPalpha IIb, GPIIb	
Immunogen	Murine platelets	
Entrez Gene ID	16399	
Gene name	Itga2b	
NCBI Full Gene Name	integrin subunit alpha 2b	
UniProt ID	Q9QUM0	
Concentration	1 mg/ml	
Preparation	Purified by protein-G affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.