

Technical Data Sheet

Product

Anti-CRP Biotin

Cat. Number/Size

1B-484-C100

0.1 mg

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CRP
Clone	C7
Format	Biotin
Reactivity	Human
Application	ELISA
Isotype	Mouse IgG1
Specificity	The antibody C7 recognizes human CRP, an 117 kDa acute phase protein; it recognizes antigen both in the presence and in the absence (samples containing EDTA) of Ca ²⁺ .
Other names	C-Reactive Protein
Immunogen	Purified human CRP

Concentration	1 mg/ml
Preparation	The purified antibody is conjugated with biotin-LC-NHS under optimum conditions. The reagent is free of unconjugated biotin.
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on vial label.

Images and References www.exbio.cz

Unless indicated otherwise, all products are For Research Use Only and not for diagnostic use. In vivo diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s. © 2019 EXBIO Praha, a.s. All rights reserved.