

Technical Data Sheet

Product	Anti-LARGE1 Purified	
Cat. Number/Size	11-702-C025	0.025 mg
	11-702-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	LARGE1
Clone	LARGE-02
Format	Purified
Reactivity	Human
Application	FC (QC tested), WB
Application details	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Isotype	Mouse IgG2b
Specificity	The mouse monoclonal antibody LARGE-02 recognizes human LARGE1, a glycosyltransferase expressed mainly in the Golgi apparatus. Crossreactivity with LARGE2 was not determined.
Other names	Xylosyl- and glucosyl transferase 1
Immunogen	Recombinant fragment of human LARGE1 (amino acids 35-142)
Entrez Gene ID	9215
Gene name	LARGE1
NCBI Full Gene Name	LARGE xylosyl- and glucuronyltransferase 1
UniProt ID	O95461

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.