

Technical Data Sheet

Product	Anti-LAT Purified	
Cat. Number/Size	11-360-C025	0.025 mg
	11-360-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	LAT
Clone	LAT-01
Format	Purified
Reactivity	Human
Negative species	Mouse
Application	WB, IP, FC
Isotype	Mouse IgG1
Specificity	The antibody LAT-01 reacts with an intracellular epitope of LAT, a 36-38 kDa transmembrane adaptor protein expressed by T cells, pre-B cells, NK cells, mast cells and platelets.
Other names	IMD52
Immunogen	Bacterially produced recombinant polypeptide corresponding to the entire cytoplasmic domain of human LAT.
Entrez Gene ID	27040
Gene name	LAT
NCBI Full Gene Name	linker for activation of T cells
UniProt ID	O43561

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.