

Technical Data Sheet

| | | |
|------------------|---|---------------|
| Product | Anti-N-sulfated heparan sulfate Purified | |
| Cat. Number/Size | 11-821-C100 | 0.1 mg |
| | For Research Use Only. | |
| | Not for use in diagnostic or therapeutic procedures. | |

| | |
|---------------------|---|
| Antigen | N-sulfated heparan sulfate |
| Clone | HepSS-1 |
| Format | Purified |
| Reactivity | Species independent |
| Application | WB, ICC, FC (QC tested), IHC, IP |
| Application details | Flow cytometry: Recommended dilution: 3-12 µg/ml. |
| Isotype | Mouse IgM |
| Specificity | The mouse monoclonal antibody HepSS-1 (also known as HepSS1) recognizes N-sulfated heparan sulfate (extracellular antigen) present in many species. |
| Other names | N-sulfo-rich heparan sulfate |
| Immunogen | MethA murine fibrosarcoma |

| | |
|----------------------|--|
| Concentration | 1 mg/ml |
| Preparation | Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods). |
| Formulation | Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide |
| Storage and handling | Store at 2-8°C. Do not freeze. |

| | |
|-----------------------|--|
| Images and References | www.exbio.cz |
|-----------------------|--|

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.