



	Technical Data Sheet	
Product	Anti-HDAC6 Purified	
Cat. Number/Size	11-842-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	HDAC6	
Clone	178	
Format	Purified	
Reactivity	Human	
Application	ICC, ELISA, IP	
Application details	Immunocytochemistry: Recommended dilution: 5-10 µg/ml.	
Isotype	Mouse IgG1 kappa	
Specificity	The mouse monoclonal antibody 178 recognizes zinc finger domain of human histone deacetylase 6 (HDAC6; an intracellular antigen). Crossreactivity with other species was not tested.	
Other names	Histone deacetylase 6	
Immunogen	human HDAC6	
Entrez Gene ID	10013	
Gene name	HDAC6	
NCBI Full Gene Name	histone deacetylase 6	
UniProt ID	Q9NSW6	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.