



	Technical Data Sheet	
Product	Anti-Hu CD11a Purified	
Cat. Number/Size	11-380-C025	0.025 mg
	11-380-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD11a	
Clone	MEM-83	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested), IP	
Application details	Flow cytometry: Recommended dilution: 1 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody MEM-83 reacts with an extracellular epitope of CD11a (alpha subunit of human LFA-1), a 170-180 kDa type I transmembrane glycoprotein expressed on B and T lymphocytes, monocytes, macrophages, neutrophils, basophils and eosinophils.	
Other names	LFA-1, LFA1A, ITGAL	
Workshop	HLDA IV: WS Code N 211	
Immunogen	Human peripheral blood lymphocytes	
Entrez Gene ID	3683	
Gene name	ITGAL	
NCBI Full Gene Name	integrin subunit alpha L	
UniProt ID	P20701	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.