



	Technical Data Sheet	
Product	Anti-Hu CD355 Purified	
Cat. Number/Size	11-910-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD355	
Clone	Cr24.1	
Format	Purified	
Reactivity	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 3-6 µg/ml.	
Isotype	Mouse IgG2b	
Specificity	The mouse monoclonal antibody Cr24.1 recognizes an extracellular epitope of human CD355, a type I transmembrane glycoprotein, expressed on subsets of CD8+ and CD4+ T cells, on activated NK cells and NKT cells.	
Other names	CRTAM	
Immunogen	CD355-P815 cells	
Entrez Gene ID	56253	
Gene name	CRTAM	
NCBI Full Gene Name	cytotoxic and regulatory T cell molecule	
UniProt ID	O95727	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.