

Technical Data Sheet

Product

Anti-Hu CD303 PE

Cat. Number/Size

1P-943-T100

100 tests

For Research Use Only.

Not for use in diagnostic or therapeutic procedures.

Antigen	CD303
Clone	15E3
Format	PE
Reactivity	Human
Application	FC (QC tested)
Application details	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Excitation laser	blue (488 nm)
Isotype	Mouse IgG1 kappa
Specificity	The mouse monoclonal antibody 15E3 recognizes an extracellular epitope of human CD303, a transmembrane glycoprotein expressed on plasmacytoid dendritic cells.
Other names	BDCA2, CLEC4C, DLEC, HECL, CLECSF7, CLECSF11, PRO34150
Immunogen	CD303 ectodomain fused with human IgG Fc domain
Entrez Gene ID	170482
Gene name	CLEC4C
NCBI Full Gene Name	C-type lectin domain family 4 member C
UniProt ID	Q8WTT0

Preparation Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Formulation Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage and handling Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.