



	Technical Data Sheet	
Product	Anti-MAP2 Purified	
Cat. Number/Size	11-253-C025	0.025 mg
	11-253-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	MAP2	
Clone	MT-01	
Format	Purified	
Reactivity	Cow, Cat, Mouse, Pig, Human	
Application	IP, WB, ICC	
Isotype	Mouse IgG1	
Specificity	The antibody MT-01 reacts with Microtubule Associated Protein 2ab (MAP2ab), an intracellular antigen.	
Other names	Microtubule-associated protein 2, MAP-2	
Immunogen	Termostabile fraction of bovine brain microtubules (whole MAP2ab molecule).	
Entrez Gene ID	4133	
Gene name	MAP2	
NCBI Full Gene Name	microtubule associated protein 2	
UniProt ID	P11137	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.