



	Technical Data Sheet	
Product	Anti-Ms CD79b FITC	
Cat. Number/Size	1F-584-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD79b	
Clone	HM79	
Format	FITC	
Reactivity	Mouse	
Negative species	Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-2 µg/ml.	
Excitation laser	blue (488 nm)	
Isotype	Armenian Hamster IgG2	
Specificity	The Armenian hamster monoclonal antibody HM79 recognizes an extracellular epitope of mouse CD79b (CD79 beta, Ig beta), a component of B cell receptor (BCR) complex.	
Other names	BCR beta, Ig-beta, B29, IGB	
Immunogen	Purified CD79a/b (alpha/beta) dimers from WEHI231 cells	
Entrez Gene ID	15985	
Gene name	CD79B	
NCBI Full Gene Name	CD79b molecule	
UniProt ID	P15530	
Concentration	0.5 mg/ml	
Preparation	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.