

Technical Data Sheet

Product	Anti-Ms Blnk Purified	
Cat. Number/Size	11-411-C025	0.025 mg
	11-411-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	Blnk
Clone	PAb (411)
Format	Purified
Reactivity	Mouse
Application	WB
Application details	Western blotting: Reducing conditions.
Isotype	Rabbit polyclonal
Specificity	The polyclonal antibody reacts with mouse SLP65 / BLNK, a cytosolic adaptor protein identified as two phosphoproteins migrating at 68 and 70 kDa in SDS/PAGE (alternatively spliced forms of human SLP65).
Other names	SLP65, BASH, Ly57
Immunogen	A fusion protein representing amino acids 171-356 of human BLNK.
Entrez Gene ID	29760
Gene name	BLNK
NCBI Full Gene Name	B cell linkerprovided
UniProt ID	Q8WV28

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.