



	Technical Data Sheet	
Product	Anti-Ms CD90.2 PE	
Cat. Number/Size	1P-443-C025	0.025 mg
	1P-443-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	CD90.2	
Clone	1aG4/C5	
Format	PE	
Reactivity	Mouse	
Negative species	Rat, Human	
Application	FC (QC tested)	
Application details	Flow cytometry: Recommended dilution: 1-2 µg/ml.	
Excitation laser	blue (488 nm)	
Isotype	Mouse IgG3	
Specificity	The mouse monoclonal antibody 1aG4/C5 recognizes alloantigen Thy-1.2 of mouse Thy-1 (CD90); it does not cross-react with Thy-1.1 alloantigen. Thy-1 is an 18-35 kDa cell surface GPI anchored glycoprotein (immunoglobulin superfamily); in mouse it is abundant on thymocytes and peripheral T lymphocytes.	
Other names	Thy-1.2	
Immunogen	thymocytes from C3H/Di mice (Thy-1.2 positive)	
Entrez Gene ID	21838	
Gene name	THY1.2	
NCBI Full Gene Name	Thy-1.2 cell surface antigen	
UniProt ID	P01831	
Concentration	0.1 mg/ml	
Preparation	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.	
Formulation	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.