



	Technical Data Sheet	
Product	Anti-Hu IL-2 Purified	
Cat. Number/Size	11-936-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	IL-2	
Clone	35C3	
Format	Purified	
Reactivity	Human	
Application	IP, FC (QC tested), ELISA, ICC, WB	
Application details	Flow cytometry: Recommended dilution: 0.5-4 µg/ml. Intracellular staining.	
Isotype	Mouse IgG2b kappa	
Specificity	The mouse monoclonal antibody 35C3 recognizes human interleukin 2 (IL-2; secreted or intracellular).	
Other names	Interleukin 2	
Immunogen	recombinant human IL-2	
Entrez Gene ID	3558	
Gene name	IL2	
NCBI Full Gene Name	interleukin 2	
UniProt ID	P60568	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.