

Technical Data Sheet

Product

**Anti-hGH Purified**

Cat. Number/Size

**11-126-C025**

**0.025 mg**

**11-126-C100**

**0.1 mg**

**For Research Use Only.**

**Not for use in diagnostic or therapeutic procedures.**

---

Antigen	hGH
Clone	GH-45
Format	Purified
Reactivity	Human
Application	IHC(P), ICC, ELISA
Isotype	Mouse IgG1
Specificity	The mouse monoclonal antibody GH-45 recognizes human growth hormone (hGH), a polypeptide hormone synthesized by acidophilic or somatotropic cells of the anterior pituitary gland. The GH-45 antibody reacts with affinity constant $3.8 \times 10^{10}$ l/mol; it does not bind human prolactin or any other pituitary hormones.
Other names	Somatotropin, Growth hormone
Immunogen	Human growth hormone.
Entrez Gene ID	2688
Gene name	GH1
NCBI Full Gene Name	growth hormone 1
UniProt ID	P01241

---

Concentration	0.5 mg/ml
Preparation	Purified by protein-A affinity chromatography
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.

---

Images and References [www.exbio.cz](http://www.exbio.cz)

---

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at [www.exbio.cz](http://www.exbio.cz). EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.