



	Technical Data Sheet	
Product	Anti-HLA-G Purified	
Cat. Number/Size	11-394-C025	0.025 mg
	11-394-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	HLA-G	
Clone	MEM-G/4	
Format	Purified	
Reactivity	Human	
Application	WB	
Application details	Western blotting: Recommended dilution: 1-2 µg/ml.	
Isotype	Mouse IgG1	
Specificity	The antibody MEM-G/4 reacts with an extracellular epitope of denaturated HLA-G heavy chain. The recognized epitope is present only on HLA-G1, -G2 and -G5 molecules. HLA-G belongs to the MHC Class I molecules (MHC Class Ib; nonclassical) and it is expressed on the surface of trophoblast cells.	
Immunogen	Recombinant human HLA-G denaturated heavy chain.	
Entrez Gene ID	3135	
Gene name	HLA-G	
NCBI Full Gene Name	major histocompatibility complex, class I, G	
UniProt ID	P17693	
Concentration	1 mg/ml	
Preparation	Purified by protein-A affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.