

Technical Data Sheet

Product	Anti-SLP76 Purified	
Cat. Number/Size	11-412-C025	0.025 mg
	11-412-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	SLP76	
Clone	PAb (412)	
Format	Purified	
Reactivity	Human	
Application	WB, IHC(P)	
Isotype	Rabbit polyclonal	
Specificity	The polyclonal antibody reacts with SLP76, a 76kDa cytosolic adaptor protein that is involved in signaling of various hematopoietic cells, such as T cells, mast cells or neutrophils; in B cells, however, it is replaced by SLP65.	
Other names	SLP-76, LCP2	
Immunogen	Bacterially expressed fusion protein representing amino acids 216-434 of human SLP76 with histidine tag	
Entrez Gene ID	3937	
Gene name	LCP2	
NCBI Full Gene Name	lymphocyte cytosolic protein 2	
UniProt ID	Q13094	
Concentration	1 mg/ml	
Preparation	Purified from rabbit serum by precipitation methods	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze. Do not use after expiration date stamped on the label.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.