

Technical Data Sheet

Product	Anti-Hu CD266 Purified	
Cat. Number/Size	11-755-C025	0.025 mg
	11-755-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	

Antigen	CD266
Clone	ITEM-4
Format	Purified
Reactivity	Mouse, Human
Application	FC (QC tested), WB, IHC(F)
Application details	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG2b
Specificity	The mouse monoclonal antibody ITEM-4 recognizes an extracellular epitope of CD266 / TWEAK R, a TNFR superfamily receptor for CD255 / TWEAK, a TNF-like weak inducer of apoptosis.
Other names	TWEAKR, FN14, TNFRSF12A
Immunogen	human CD266-transfected P815 cells
Entrez Gene ID	51330
Gene name	TNFRSF12A
NCBI Full Gene Name	TNF receptor superfamily member 12A
UniProt ID	Q9NP84

Concentration	1 mg/ml
Preparation	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage and handling	Store at 2-8°C. Do not freeze.

Images and References www.exbio.cz

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.