



Technical Data Sheet

Product	Anti-Ag85b (M. tuberculosis) Purified	
Cat. Number/Size	11-189-C100	0.1 mg
	For Research Use Only.	
	Not for use in diagnostic or therapeutic procedures.	
Antigen	Ag85b (M. tuberculosis)	
Clone	PAb (189)	
Format	Purified	
Reactivity	Mycobacterium tuberculosis	
Application	WB, ELISA	
Isotype	Rabbit polyclonal	
Specificity	The rabbit polyclonal antibody to Ag85b reacts with the Mycobacterium tuberculosis protein Ag85b (FbpB), which is a diacylglycerol acyltransferase/mycolyltransferase.	
Other names	FbpB, mycolyl transferase 85B	
Immunogen	Recombinant Mycobacterium tuberculosis Ag85b produced in E. coli	
Entrez Gene ID	885785	
Gene name	fbpB	
NCBI Full Gene Name	diacylglycerol acyltransferase/mycolyltransferase	
UniProt ID	I6YBV4	
Concentration	1 mg/ml	
Preparation	Purified by ligand affinity chromatography.	
Formulation	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide	
Storage and handling	Store at 2-8°C. Do not freeze.	
Images and References	www.exbio.cz	

The product is intended For Research Use Only. Diagnostic or therapeutic applications are strictly forbidden. Products shall not be used for resale or transfer to third parties either as a stand-alone product or as a manufacture component of another product without written consent of EXBIO Praha, a.s. EXBIO Praha, a.s. will not be held responsible for patent infringement or any other violations of intellectual property rights that may occur with the use of the products. Orders for all products are accepted subject to the Term and Conditions available at www.exbio.cz. EXBIO, EXBIO Logo, and all other trademarks are property of EXBIO Praha, a.s.